



**CORFLEX**

**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



## ARM SLINGS

### TRICOT ARM SLING

- > Features thumb loop to prevent wrist drop and Velcro closure at top of sling envelope to provide a secure fit
- > Indicated for injuries to the arm, shoulder, wrist or hand and cast support to maintain arm in neutral position
- > Manufactured from soft yet durable tricot material



#### PRODUCT NUMBER

XS	S	M	L	XL
21-5060	21-5061	21-5062	21-5063	21-5064

#### SIZING GUIDE

	XS	S	M	L	XL
Width/Height	12" x 5"	15" x 6"	17" x 7"	19" x 8"	21" x 9"

Measurements are sling pouch width x height

### APPLICATION INSTRUCTIONS

1. Open sling and place arm in inside, making sure elbow is in sling as far back as possible.
2. Thread strap through d-rings.
3. Secure hook and loop so elbow is positioned at 90-degree angle. (If necessary, there is a one-time hook and loop adjustment behind the elbow to lengthen or shorten strap length as needed.)
4. Position elbow (shoulder) strap and pull behind your back.
5. Thread strap through buckle and pull tight to prevent arm from moving away from body.

### CARE & CLEANING

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

### HOW TO PLACE AN ORDER

The Corflex Customer Service Team is ready to assist you in placing Your order. You may place your orders via a variety of methods:



**TELEPHONE**

1-800-426-7353, Ext1



**FAX**

1-603-623-4111



**E-MAIL**

Service@corflex.com



**MAIL**

CORFLEX, Inc.  
Attn: Customer Service  
669 East Industrial Park Dr  
Manchester, NH 03109 USA

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - www.corflex.com

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.