



**CORFLEX**

**CORFLEX IS A GLOBALLY BASED  
MANUFACTURER AND DISTRIBUTOR OF QUALITY  
ORTHOPEDIC AND REHABILITATION PRODUCTS**



**KNEE MISC**

**CRYO PNEUMATIC KNEE SPLINT W/STAYS**

- > Anatomical knee wrap designed to combine the benefits of pneumatic compression and rigid immobilization
- > Indicated for post-op immobilization, acute injury and as part of a physical therapy or athletic training program; helps diminish swelling / edema, reduce muscle spasms and decrease pain
- > Manufactured from plush foam laminate that is Velcro compatible with rigid aluminum stays



**PRODUCT NUMBER**

<b>0 GELS</b>	<b>1 GEL</b>	<b>2 GELS</b>
<b>72-7170</b>	<b>72-7171</b>	<b>72-7172</b>

**SIZING GUIDE**

Universal  
Additional gels sold separately:  
72-7179 - Single Knee Gel

**APPLICATION INSTRUCTIONS**

[Click here.](#)

**CARE & CLEANING**

Corflex manufactures a quality line of products and to help your care for the product properly we provided the following instructions for care & cleaning:

Close all contact closures and hand wash in warm water with mild soap. Product should be air-dried, do not use high temperatures to wash or dry.

669 East Industrial Park Drive - Manchester, NH 03109  
Phone (603) 623-3344 - Fax (603) 623-4111 - [www.corflex.com](http://www.corflex.com)

© Corflex | The advice and information contained in this document is not intended as a substitute for medical counseling. If you have any questions about treatment for a specific condition, please check with a health care professional.